



## CANOLA SEED ANALYZER by MAGNETIC RESONANCE



### Simultaneous and Automatic Determination of **MOISTURE** and **TOTAL OIL** Content



- **Sample Size: 45 cm3 Bulk or Single Seed.**
- **Non-destructive Testing.**
- **No Sample Preparation.**
- **Low Maintenance.**
- **Simple Graphical Interface.**
- **Portable and flexible device for using in other applications.**

	Accuracy	Time (secs)
<b>Total Oil %</b>	<b>± 0,1%</b>	<b>4</b>
<b>Moisture %</b>	<b>±0,1%</b>	<b>4</b>

### SLK NMR Seed v3.0 Software



### Main Advantages

- Results on screen in a few seconds.
- Simple graphical interface, secure data acquisition and simultaneous processing several components.
- Automatic determination of the sample weight, reducing operating mistakes.
- It allows the import and export of data, reports generation and measurements record.
- It allows to renew calibrations, to add or to deactivate points, to show the adjustments and statistical parameters. Graphic visualization of the results and visualization of residues.
- Report for feat technical support and on-line diagnosis of the spectrometer state.
- Scanning of measurements condition and temperature.
- It allows the users management, permission of accesses, permission to modify information, etc.